



SPME Fiber Holder

- ◆ For Manual Sampling
- ◆ For Varian/CTC CombiPal/Gerstel MPA 2/Thermo TriPlus Autosampler
- ◆ SPME Portable Field Sampler



Application for SPME

- ◆ Environmental analysis of water samples
- ◆ Headspace analysis of trace impurities solid ppt odor analysis
- ◆ Flavor analysis of food products
- ◆ Forensic analysis of arson/explosives samples
- ◆ Toxicology analysis for blood alcohol or drugs in urine/serum
- ◆ Surfactants, other industrial applications

SPME Accessories

- ◆ SPME sampling stand
- ◆ SPME inlet guide
- ◆ SPME inlet liner
- ◆ Heater block
- ◆ Vials for SPME sampling stand

SPME Fiber Assemblies

For manual holder/autosampler/ Multi Fiber Exchanger

- ◆ Polydimethylsiloxane (PDMS)
- ◆ Carboxen/Polydimethylsiloxane (CAR/PDMS)
- ◆ Polydimethylsiloxane/Divinylbenzene (PDMS/DVB)
- ◆ Polyacrylate (PA)
- ◆ Divinylbenzene/Carboxane/Polydimethylsiloxane (DVB/CAR/PDMS)
- ◆ Carbowax-Polyethylene Glycol (PEG) coating

