

Automated Drying & Concentration System





DryVap Concentrator System



- ◆ Concentrate up to six (6) sample extracts simultaneously or individually
- ◆ Concentrate using a combination of heat, vacuum and nitrogen sparge gas
- ◆ Automatic endpoint detection and shut-off
- ◆ Automatic in-line removal of residue water using DryDisk™ technology
- ◆ Sealed evaporation tubes to prevent loss of analytes
- ◆ Superior analyte recovery due to Kuderna-Danish (KD) effect

DryDisk™ Separation Membrane

- ◆ Removal of residual water from extracting solvent using Gore-Tex Technology
- ◆ Eliminate the need of a drying agent



200mL Evaporation Tube

- ◆ Graduated tip in various volume (0.4, 0.9, 1.0 and 2.0 mL)
- ◆ Direct to GC Vial
- ◆ Solid tip for solvent removal

