


The World's Pre-eminent Manufacturer of Thermal Desorption Tubes  
Trusted by US Military Since 1991



## Thermal Desorption Tubes

- ◆ Inert Stainless Steel
- ◆ Stainless Steel (Surgical grade 316L)
- ◆ Glass
- ◆ PTFE

Standard and custom-made Thermal Desorption tubes for ATD systems from:

- ◆ PerkinElmer
- ◆ CDS Dynatherm
- ◆ Gerstel
- ◆ Teledyne Tekmar
- ◆ Markes
- ◆ Shimadzu
- ◆ Dani
- ◆ And more ...

## COGNI™ System

Every Camsco thermal desorption tube is equipped with COGNI™ System

- ◆ **Serial numbers** and **barcodes** for storage, inventory management and data retrieval
- ◆ **Circumferential bands** for quick sorbent identification
- ◆ Complete family of **micro caps** to maintain and ensure the integrity of samples



## Application

Camsco provides different types of tubes for specific application:

- ◆ Tenax TA Tube
- ◆ Air Toxics Tube
- ◆ Universal Tube
- ◆ TO -17 Tube Style 1, Style 2 & Style 3
- ◆ Cold Traps
- ◆ Benzene Monitoring Tube
- ◆ BTEX Tube
- ◆ BTEX Capture System
- ◆ NIOSH 2549 Tube
- ◆ Soil Gas Tube
- ◆ Soil Vapor Intrusion SVI Tube
- ◆ Universal Tube